



Know VTE

APPLICATION GUIDE



AGENDA

- Introduction
- How to install?
- Know VTE application guide (How to use?)
- Q&A



INTRODUCTION




HOW TO INSTALL?

1. Download from the provided links
2. On iOS, Trust the developer license
3. Use the application

1- DOWNLOAD

- Open the following link on your cell phone or scan the QR code
 - [Download Link](#)
 - Choose your platform by clicking on the icon (iOS or Android)





INSTALLATION ON IOS

- 1- After installation, search for the "Know VTE" app icon.
- 2- When you try to run the app for the first time you will probably get the "Untrusted Enterprise Developer"
- 3- Click cancel and open your device settings
- 4- Go to General
- 5- Scroll until you find the "Device management" option and click on it
- 6- Click on "Eye Egypt"
- 7- Then click on Trust "Eye Egypt"
- 8- Finally click on Trust
- 9- Now you can run the app :)

INSTALLATION ON IOS

1- After installation, search for the "Know VTE" app icon



2- When you try to run it for the first time, you will probably get the following message:



3- Click cancel and go to settings




INSTALLATION ON IOS


4- Go to General




5- Scroll until you find Device management and click it



6- Click on Eye Egypt




7- Then click on Trust "Eye Egypt"



INSTALLATION ON IOS

8- Finally click on **Trust**




9- Now you can open the application normally :)



KNOW VTE APPLICATION GUIDE (HOW TO USE?)

IN THIS GUIDE WE WILL PRESENT THE VTE MOBILE APP FEATURES AND ELABORATE ON HOW IT WORKS.



SIGN IN

Email
dr.test@test-hos-one.com


Password

Remember Me

SIGN IN →

Don't have an account? [Sign Up!](#)

sanofi |  | 



SIGN UP

Welcome

Email *

Full Name *

Phone Number *

Password *


Confirm Password *

SIGN UP AND LOGIN

- Predefined Users (doctors accounts added by the admin) can login to the application using their credentials found in the received "VTE - Registration" email.
- Other doctors can sign up to use the application's features by:
 - Going to the sign-up page
 - Fill in the Email address, full name and password then click on sign up, you will receive a confirmation email which contains a button that when clicked your account will be verified and ready to use.


FIND/ADD PATIENT

- After logging in, one of the first actions you can do is to find an old registered patient or add a new one.
- Using the upper search field, you can search for the patient with his/her national id or name.
- Or you can add a new patient by clicking on the "Add patient" button, you will be redirected to the "Create patient profile" page.



FIND/ADD PATIENT

- In the create patient profile form, the doctor can add the following fields:
 - الإسم الرباعى باللغة العربية
 - Full name in English
 - Gender
 - Age Range
 - Medical Record ID
 - National ID

A mobile application form titled 'Create Patient Profile' with a blue and red header. The form is titled 'ADD PATIENT' and contains several input fields and buttons. The fields are: 'الإسم بالعربية (رباعى) *' (Arabic name, 4 characters), 'Full Name in English', 'Gender *' (with 'Male' and 'Female' buttons), 'Age Range *' (with buttons for '18-30', '31-40', '41-50', '51-60', and '61+'), and 'Medical Record ID'.

Create Patient Profile

ADD PATIENT

* الإسم بالعربية (رباعى)

Full Name in English

Gender *

Male Female

Age Range *


18-30 31-40 41-50

51-60 61+

Medical Record ID


PATIENT PROFILE

- Once the profile is created you will automatically get redirected to the patient's profile where you can start adding assessment(s)
- You can also go to the patient profile from the homepage by typing the patient's name in the search field.
- The patient profile contains the following info:
 - Patient's name in Arabic and English languages
 - Gender
 - Medical record ID
 - List of past assessments (if the patient hadn't had any assessments, it will appear empty)

A mobile app interface for a patient profile. The header is a blue-to-red gradient bar with a hamburger menu icon and the text 'Patient Profile'. Below the header, the patient's name is displayed in Arabic: 'محمد أحمد حسن سيد'. There are two light blue boxes: one for 'Patient Name En' containing 'Mohamed Ahmed Hassan Said', and one for 'Gender' containing 'Male'. Below these is a box for 'Medical Record ID' containing '9911552'. Underneath is the section 'Assessment List' which is currently empty. At the bottom, there is a blue button labeled 'ADD ASSESSMENT' with two radio button options: 'First Assessment' and 'Reassessment'. At the very bottom is a red button labeled 'Start' with a right-pointing arrow.


ADD ASSESSMENT

- In the patient profile, you can add an assessment by:
 - First choose whether it's a first assessment or reassessment.
 - If it's a reassessment you must choose a reason
- For the first time, after creating the patient profile, you can create a first assessment only.
- Select first assessment or reassessment then click "Start".

A screenshot of a mobile application interface for a patient profile. The title bar is blue and red with a white hamburger menu icon and the text 'Patient Profile'. Below the title bar, the patient's name is displayed in Arabic: 'محمد أحمد حسن سيد'. There are two input fields: 'Patient Name En' containing 'Mohamed Ahmed Hassan Said' and 'Gender' containing 'Male'. Below these is a 'Medical Record ID' field with the value '9911552'. Underneath is a section titled 'Assessment List' which is currently empty. At the bottom, there is a section titled 'ADD ASSESSMENT' with two radio button options: 'First Assessment' and 'Reassessment'. At the very bottom is a red button with the text 'Start' and a right-pointing arrow.


PICK ASSESSMENT TYPE

- Next, you will choose a specialty which will determine the type of assessment
- The Gynecology assessment will not be available for male patients.
- The assessments are as the following:
 - Oncology (Khorana and HasBled assessments)
 - Medical (Caprini assessment)
 - Surgical (Caprini assessment)
 - Gynecology (RCOG assessment)




ASSESSMENT

- Then you should answer all the questions in the assessment.
- Each question will have multiple answers, and it also might contain a small description at the bottom to elaborate more about the question meaning an objective.

A screenshot of a mobile application interface for a "Khorona Assessment". The screen has a dark blue header with a white back arrow and the text "Khorona Assessment". Below the header, the main content area is white and contains the following text: "Know the patient's Risk for Blood Clots & Save his life", "Answer These Simple Questions", and a progress indicator "1/5" with a red progress bar. The question is "Cancer type" and lists several options with radio buttons: Stomach, Pancreas (selected), Lung, Lymphoma, Gynecologic, Bladder, Testicular, and Other. At the bottom of the screen, there are two navigation buttons: a white back arrow and a red forward arrow.


GET ASSESSMENT RESULTS

- After answering all the questions, you will find the assessment score and the recommended prophylaxis
- You must then choose whether you are going to follow the guidelines, or try a different approach
- If you chose to follow a different approach, you should mention the justification



GET ASSESSMENT RESULTS

- After submitting the assessment, you will get a finishing message with two buttons.
- You can either finish the assessment or show the full report.
- If you choose to show the full report, and for the first time only, you will get a message asking you to allow access to media files, so the Know VTE application can save the report on your device.
- Click "Allow", then the report will appear.



FULL REPORT

- The full assessment report will include:
 - The date and time of the assessment
 - Doctor name
 - Patient's name, gender, and medical record ID
 - All the questions of the assessment with their answers and scores
 - Final assessment scores
 - Recommended prophylaxis
 - Variance
 - Justification comment (If exist)

Assessment Report
✕



المملكة العربية السعودية
وزارة الصحة والرفاهية
KSAHES

Date:
Sep 11, 2022

Time:
5:38 PM

Doctor Name: Mohamed Medhat

Medical Record ID:
9911552

Patient Name: محمد أحمد حسن سيد

Gender: Male

Khorona Questions	Answer	Score
Cancer type	Gynecologic	1
Pre-chemotherapy platelet count $\geq 350 \times 10^9/L$	Yes	1
Hemoglobin level < 10 g/dL or using RBC growth factors	Yes	1
Pre-chemotherapy leukocyte count $> 11 \times 10^9/L$	Yes	1
BMI ≥ 35 kg/m ²	Yes	1

HasBled Questions

AnswerScore


Save

Print

Share


FULL REPORT

- You will have the ability to Save, Print, or Share the report.
- If you click on Save, the report will be saved on your mobile storage and another success message will appear with an open the saved file button.
- Click on "OPEN FILE" to view the full report in PDF format.
- A PDF reader must be installed on your device to be able to view the full report, you can download the PDF reader from here:
 - For android devices: [Download](#)
 - For iPhone devices: [Download](#)




FULL REPORT

- If you clicked on Print, you would get prompted that a printer with Wi-Fi connection is needed in order to print directly from your mobile.
- If you have a printer connected to your Wi-Fi network, then click on "I Understand".
- After that you can, choose your printer and click "Print" to print the full report.




FULL REPORT

- You will also be able to share the report using the "Share" button.
- Just click on the "share" button, and you can share the report in PDF format through WhatsApp, email, Messenger, or any other application that you prefer.



FULL REPORT

- You can close the report by clicking on the “x” sign at the top right part of the page.
- You will be redirected back to the patient’s profile where you can find all the assessments of the patient and you can open the report of any of the assessments by clicking on “Open full report button”.




HOW TO USE? (IN NUTSHELL)




CONTACT US

- You can contact us regarding any issue you face during your experience with the “Know VTE” app through the contact us page.
- You can find the contact us page from the top left menu 
- Click on “Contact Us”, choose the topic and enter your message, then click on “Submit”.
- An email will be sent to Gahar with your message and Gahar’s representative will contact you as soon as possible.



SETTINGS AND NOTIFICATIONS

- Finally, you can visit the settings from the side menu, where you can disable/enable receiving notifications and Dark mode
- The notifications inbox contains a list of received notifications





THANK YOU